
# ***BIARTICULAR***

***Cabeza Biarticular***  
***Biarticular Head***

 **surgival**  
GRUPO COSÍAS

# Cabeza Biarticular

## Biarticular Head

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The **Biarticular Head**, also named **Bipolar or Bicentric**, was conceived fundamentally to reduce the friction that developed, in most cases and in an advanced form, in the cartilage acetabular to the direct contact of the metallic head of an endoprosthesis in a classic hemi arthroplasty.

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The **Biarticular Head** consists of two parts: a metallic **Cup** that is in contact with the acetabulum and one **Insert** in polyethylene that adapts inside the metallic cup whose interior fixes the head of the femoral stem. In this way there are two articulations made, one between acetabulum and cup, and another between the inset and femoral head.

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The **Cup** is made of stainless steel (AISI 316 LVM), and is presented in 15 sizes from 41 to 55 mm, with a cadence of 1 mm, in function of the diameter of the femoral head.

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The **Insert** is made of polyethylene of a ultra high molecular weight (UHMWPE), that is translated in a minor "debri" by stress. It is presented in four sizes with an internal diameter of 28 mm.

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



1 · The way of beginning is to the discretion of the surgeon.

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4 · Opening and exploration of the medular cavity with a narrow and long spoon.

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5 · Introduction of the rasp in order to make the femoral cavity of the femoral prosthesis.

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7 · On the definitive stem we can carry out a trial reduction using the trial cups and trial heads, the instrumentation set includes trial 28 mm diameter heads and taper 12/14.

**Note:** If the chosen stem is SURGIVAL, we will be able to carry out previously trial reductions on heads and cups even on the rasp.

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8 · Assembling the different components following the next steps:


1<sup>st</sup> Step: To introduce the head of 28 mm diameter on pressure in the polyethylene insert.

2<sup>nd</sup> Step: The block head-insert is introduced in the Biarticular cup, check it is correctly fit together.

**Note:** The assembling of the femoral head, insert and cup is advisable to carry out on the instrumentation table. In the cases that are considered opportune it is also possible to assemble directly on the neck of the femoral stem.

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9 · Impactation of the assembling Biarticular components on the stem and articulation reduction.

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# Instrumental Instruments

Set completo instrumental biarticular  
Biarticular instrumentation complete set Ref. A1519000

Bandeja instrumental  
Instrumentation tray Ref. A1519011

1	Cabeza prueba biarticular Trial biarticular head
Ref. A1519541	Ø 41 mm
Ref. A1519542	Ø 42 mm
Ref. A1519543	Ø 43 mm
Ref. A1519544	Ø 44 mm
Ref. A1519545	Ø 45 mm
Ref. A1519546	Ø 46 mm
Ref. A1519547	Ø 47 mm
Ref. A1519548	Ø 48 mm
Ref. A1519549	Ø 49 mm
Ref. A1519550	Ø 50 mm
Ref. A1519551	Ø 51 mm
Ref. A1519552	Ø 52 mm
Ref. A1519553	Ø 53 mm
Ref. A1519554	Ø 54 mm
Ref. A1519555	Ø 55 mm

2	Cabeza prueba para vástago Ø 28 mm Trial head for stem Ø 28 mm	
	C/C cono 12/14 Short Neck 12/14 taper	Ref. A1536040
	C/M cono 12/14 Medium Neck 12/14 taper	Ref. A1536041
	C/L cono 12/14 Long Neck 12/14 taper	Ref. A1536042
	C/XL cono 12/14 Extra Long Neck 12/14 taper	Ref. A1536043

3	Discos probadores (par/impar) Tester disc (even/odd)	Ref. A1519020
4	Posicionador extractor Positioner-extractor	Ref. A1519021
5	Cinzel hueco Mortising chisel	Ref. A1160020
6	Extractor P/Insertos biarticular Extractor for biarticular insert	Ref. A1519022

Caja vacía 1 Fil 60 x 30 x 11  
Empty container 1 Fil 60 x 30 x 11 Ref. B4500001


CE0318

*Tecnología Europea de Vanguardia*  
*Advanced European Technology*

Distribuido por: **Distributed by:**

Parque Tecnológico - C/ Leonardo Da Vinci, 12-14 · 46980 Paterna · (Valencia, Spain)  
Tel. (+34) 96 131 80 50 · Fax (+34) 96 131 80 95 · E-mail: [surgival@surgival.com](mailto:surgival@surgival.com)  
[www.surgival.com](http://www.surgival.com)